IPL-Sq

# Presentation IPL-Sq®

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# CASES AND RESULTS

## HAIR REMOVAL

Photo from VYDENCE Training Center  
São Paulo, SP

IPL: Filter 640 nm, 16 J/cm<sup>2</sup>, 30 ms.

# CASES AND RESULTS

## HAIR REMOVAL

Photo from VYDENCE Training Center  
São Paulo, SP


IPL: Filter 640 nm, 17 J/cm<sup>2</sup>, 30 ms

# CASES AND RESULTS

## HAIR REMOVAL

Photo from VYDENCE Training Center  
São Paulo, SP

IPL: Filter 640 nm, 17 J/cm<sup>2</sup>, 20 ms

# CASES AND RESULTS

## HAIR REMOVAL

Photo from VYDENCE Training Center  
São Paulo, SP

**BEFORE**

**AFTER**

IPL: Filter 640 nm, 16 J/cm<sup>2</sup>, 30 ms

# CASES AND RESULTS

## HAIR REMOVAL

Photo from VYDENCE Training Center  
São Paulo, SP

**BEFORE**

**AFTER**

IPL: Filter 640 nm, 16–17 J/cm<sup>2</sup>, 30 ms, 8 sessions.

# CASES AND RESULTS

## HAIR REMOVAL

Photo from VYDENCE Training Center  
São Paulo, SP

**BEFORE**

**AFTER**

IPL: Filter 640 nm, 16–17 J/cm<sup>2</sup>, 30 ms, 8 sessions.

# CASES AND RESULTS

## HAIR REMOVAL

Photo from VYDENCE Training Center  
São Paulo, SP

**BEFORE**

**AFTER**

IPL: Filter 640 nm, 16 J/cm<sup>2</sup>, 30 ms

# CASES AND RESULTS

## HAIR REMOVAL

Photo from VYDENCE Training Center  
São Paulo, SP


IPL: Filter 640 nm, 16 J/cm<sup>2</sup>, 30 ms

# CASES AND RESULTS

## PIGMENTED LESIONS

Photo from VYDENCE Training Center  
São Paulo, SP

**BEFORE**

**AFTER**

IPL: Filter 580 and 540 nm, 15–16 J/cm<sup>2</sup>, 15 ms, 3 sessions.

# CASES AND RESULTS

## PIGMENTED LESIONS

Photo from VYDENCE Training Center  
São Paulo, SP

IPL: Filter 580 and 540 nm, 15–17 J/cm<sup>2</sup>, 15 ms, 3 sessions.

# CASES AND RESULTS

## PIGMENTED LESIONS

Photo from VYDENCE Training Center  
São Paulo, SP

IPL: Filter 580 and 540 nm, 15–16 J/cm<sup>2</sup>, 15 ms, 3 sessions.

# CASES AND RESULTS

## PIGMENTED LESIONS

Photo from VYDENCE Training Center  
São Paulo, SP

IPL: Filter 580 and 540 nm, 15–16 J/cm<sup>2</sup>, 15 ms, 3 sessions.

# CASES AND RESULTS

## PIGMENTED LESIONS

Photo from VYDENCE Training Center  
São Paulo, SP

IPL: Filter 640 nm, 14 J/cm<sup>2</sup>, 30 ms + Filter 540 nm, 15-16 J/cm<sup>2</sup>, 15 ms + AR 5%

# CASES AND RESULTS

ipl·sq

## VASCULAR LESIONS

Photo from VYDENCE Training Center  
São Paulo, SP

**BEFORE**

**AFTER**


IPL: Filter 580 and 540 nm, 20 J/cm<sup>2</sup>, 15 ms.

# CASES AND RESULTS

## VASCULAR LESIONS

Photo from VYDENCE Training Center  
São Paulo, SP

**BEFORE**

**AFTER**

IPL: Filter 580 and 540 nm, 20 J/cm<sup>2</sup>, 15 ms..

# CASES AND RESULTS

## VASCULAR LESIONS

Photo from VYDENCE Training Center  
São Paulo, SP

**BEFORE**

**AFTER**

IPL: Filter 580 and 540 nm, 20 J/cm<sup>2</sup>, 15 ms.

# CASES AND RESULTS

## ACNE VULGARIS

Photo from VYDENCE Training Center  
São Paulo, SP

IPL: Filter 640 nm, 12 J/cm<sup>2</sup>, 100 ms + Filter 400 nm, 10 J/cm<sup>2</sup>, 30 ms..

# CASES AND RESULTS

## ACNE VULGARIS

Photo from VYDENCE Training Center  
São Paulo, SP

**BEFORE**

**AFTER**

IPL: Filter 640 nm, 12 J/cm<sup>2</sup>, 100 ms + Filter 400 nm, 10 J/cm<sup>2</sup>, 30 ms..

# CASES AND RESULTS

## RECENT SCARS

Photo courtesy of Dr. Carlos Alberto Ferreira  
Campinas, SP

ANTES

DEPOIS

IPL: Filter 540 nm, 10-15 J/cm<sup>2</sup>, 10-20 ms, 3 sessions.

1. Klein A, Steinert S, Baeumler W, Landthaler M, Babilas P. PHOTOEPILATION WITH A DIODE LASER VS. INTENSE PULSED LIGHT: A RANDOMIZED, INTRAPATIENT LEFT-TO-RIGHT TRIAL. *J Dermatol*. 2013 Jun;168(6):1287-93.
2. Haak CS, Nymann P, Pedersen AT, Clausen HV, Feldt Rasmussen U, Rasmussen AK, Main K, Haedersdal M. HAIR REMOVAL IN HIRSUTE WOMEN WITH NORMAL TESTOSTERONE LEVELS: A RANDOMIZED CONTROLLED TRIAL OF LONG-PULSED DIODE LASER VS. INTENSE PULSED LIGHT. *Br J Dermatol*. 2010 Nov;163(5):1007-13.
3. Amin SP, Goldberg DJ. CLINICAL COMPARISON OF FOUR HAIR REMOVAL LASERS AND LIGHT SOURCES. *J Cosmet Laser Ther*. 2006 Jun;8(2):65-8.
4. Wu DC, Friedmann DP, Fabi SG, Goldman MP, Fitzpatrick RE. COMPARISON OF INTENSE PULSED LIGHT WITH 1,927-NM FRACTIONATED THULIUM FIBER LASER FOR THE REJUVENATION OF THE CHEST. *Dermatol Surg*. 2014 Feb;40(2):129-33.
5. Chan CS, Saedi N, Mickle C, Dover JS. COMBINED TREATMENT FOR FACIAL REJUVENATION USING AN OPTIMIZED PULSED LIGHT SOURCE FOLLOWED BY A FRACTIONAL NON-ABLATIVE LASER. *LASERs Surg Med*. 2013 Sep;45(7):405-9.
6. Ooe M, Seki T, Miura T, Takada A. COMPARATIVE EVALUATION OF WRINKLE TREATMENTS. *Aesthetic Plast Surg*. 2013 Apr;37(2):424-33.
7. Thaysen-Petersen D1, Lin JY, Nash J, Beerwerth F, Wulf HC, Philipsen PA, Haedersdal M. THE ROLE OF NATURAL AND UV-INDUCED SKIN PIGMENTATION ON LOW-FLUENCE IPL-INDUCED SIDE EFFECTS: A RANDOMIZED CONTROLLED TRIAL. *LASERs Surg Med*. 2014 Feb;46(2):104-11.
8. Srinivas CR, Kumaresan M. LASERS FOR VASCULAR LESIONS: STANDARD GUIDELINES OF CARE. *Indian J Dermatol Venereol Leprol*. 2011 May-Jun;77(3):349-68.
9. Goldberg DJ. LASER REMOVAL OF PIGMENTED AND VASCULAR LESIONS. *J Cosmet Dermatol*. 2006 Sep;5(3):204-9.
10. Liu LH, Fan X, An YX, Zhang J, Wang CM, Yang RY. RANDOMIZED TRIAL OF THREE PHOTOTHERAPY METHODS FOR THE TREATMENT OF ACNE VULGARIS IN CHINESE PATIENTS. *Photodermatol Photoimmunol Photomed*. 2013 Dec 9
9. Taylor M, Porter R, Gonzalez M. INTENSE PULSED LIGHT MAY IMPROVE INFLAMMATORY ACNE THROUGH TNF-A DOWN-REGULATION. *J Cosmet Laser Ther*. 2014 Apr;16(2):96-103
10. Raulin C, Greve B, Grema H. IPL TECHNOLOGY: A REVIEW. *LASERs Surg Med*. 2003;32(2):78-87.
11. Wat H, Wu DC, Rao J, Goldman MP. APPLICATION OF INTENSE PULSED LIGHT IN THE TREATMENT OF DERMATOLOGIC DISEASE: A SYSTEMATIC REVIEW. *Dermatol Surg*. 2014 Apr;40(4):359-77.
12. Town G, Ash C, Eadie E, Moseley H. MEASURING KEY PARAMETERS OF INTENSE PULSED LIGHT (IPL) DEVICES. *J Cosmet LASER Ther* 2007;9:148-160.
13. Clarkson DM, Gill D, Odeke M. A TIME-RESOLVED INTENSE PULSED LIGHT SPECTRAL ANALYSIS SYSTEM. *LASERs Med Sci* 2008;23:59-64; 16. Ash C, Town G, Bejorring P. RELEVANCE OF THE STRUCTURE OF TIME-RESOLVED SPECTRAL OUTPUT TO LIGHT-TISSUE INTERACTION USING INTENSE PULSED LIGHT (IPL). *LASERs in Surg and Med* 2008;40:83-92

ipl·są

 **vyidence**  
medical

 etherea™

**Thank you**

---

Marta Osóbka Panda Medical by Panda Group

proprietary and confidential