



The image shows the Etherea MX Astanza Laser device, a white and grey medical instrument. It features a large, tilted color touchscreen at the top displaying the 'etherea' logo and a colorful abstract graphic. Below the screen is a circular control panel with a red button. The device is mounted on a stand with various adjustment knobs and handles.

# Etherea's MX Astanza Laser

Dallas, October 2018

*prepared and reviewed by* Lutécia Costa, OCT 2018



A large, stylized blue graphic element consisting of several overlapping, curved shapes that resemble a funnel or a large letter 'A', positioned on the right side of the slide.

**Lutécia Ferro Costa**  
**Internacional Trade Manager**  
**[lcosta@vydence.com](mailto:lcosta@vydence.com)**

# LEAD RSHIP

#1 LASER PLATFORM IN THE WORLD'S 2ND LARGEST AESTHETIC MARKET.

**MORE SYSTEMS**  
SOLD THAN ALL  
OTHERS COMBINED\*

**750,000**  
TREATMENTS TO  
DATE AND GROWING

**93%+**  
OVERALL USER'S  
SATISFACTION


SOLD IN ALL MAJOR  
MARKETS **AROUND**  
**THE WORLD**

VERSATILITY:  
UP TO **7** DIFFERENT  
TECHNOLOGIES

**70+**  
TRATAMENT INDICATIONS  
AVAILABLE

# ETHEREA®: Key Product Milestones

CONTINUOUS IMPROVEMENT OVER THE PAST 10 YEARS



# ETHEREA®: Research and Clinical Development

4

APPROVED BY THE BEST MEDICAL INSTITUTIONS OF BRAZIL



**support to clinical research and development in PUC–Campinas (SP), UNICAMP (SP), HCFMUSP Ribeirão Preto (SP), FAMERP (SP), HCFMUSP São Paulo (SP), Santa Casa de São Paulo Einstein (SP), Santa Casa de Porto Alegre (RS), Lauro de Souza Lima Institute (SP), FMJ Hospital (SP), HRAN Brasília (DF), Hospital das Clinicas de Porto Alegre (RS), Santa Casa de Misericórdia of Rio de Janeiro (RJ), and others;**

142

### Original Articles

**Authors:**

Célia Luiza Kalil<sup>1</sup>  
Valéria Campos<sup>2</sup>  
Clarissa Prieto Herman Reinehr<sup>3</sup>  
Christine Rachelle Prescendo Chaves<sup>4</sup>

<sup>1</sup> Preceptor and Head of the Cosmiatry Ambulatory, Dermatology Service, Santa Casa de Misericórdia de Porto Alegre (RS), Brazil.

<sup>2</sup> Dermatologist physician at private practice - Jundiaí (SP), Brazil; Post graduate degree in Dermatology and Laser, Harvard Medical School - Massachusetts, USA.

<sup>3</sup> Dermatologist physician at private practice - Porto Alegre (RS), Brazil.

## *Laser toning e drug delivery: estudo-piloto utilizando laser Q-switched Nd:YAG 1064nm*

*Laser toning and drug delivery: a pilot study using laser Q-switched laser 1064nm*

DOI: <http://dx.doi.org/10.5935/scd1984-8773.201682727>

**ABSTRACT**

**Introduction:** Laser toning technique is performed with the Q-switched Nd:YAG laser 1064 nm, aiming at stimulating neocollagenesis. The technique can also be associated with the application of suitable medicines for drug delivery, increasing its potential for skin permeation.

**Objective:** To evaluate the results of laser toning Q-switched Nd:YAG laser 1064 nm associated with drug delivery through a pilot study.

**Methods:** Four patients underwent four laser sessions with application of a formula for drug delivery or placebo, fortnightly.

**Results:** According to the photographic assessment, laser toning promoted improvement of

# ETHEREA®: Key Opinion Leaders

USED AND APPROVED BY REFERRAL DOCTORS



Rodrigo Kikuchi, MD  
Brasil



Eliandre Palermo, MD  
Brasil



Sergio Talarico, MD  
Brasil



Luiza Pitassi, MD  
Brasil



Valéria Campos, MD  
Brasil



Célia Kalil, MD  
Brasil



Alexandre Filippo, MD  
Brasil



Doris Hexsel, MD  
Brasil



Emerson Lima, MD  
Brasil



Nuno Osório, MD  
Brasil



Denise Steiner, MD  
Brasil



Emerson Alves, MD  
Brasil



Maurice Adatto, MD  
Switzerland



Ricardo Shiratsu, MD  
Brasil



Elisete Crocco, MD  
Brasil



Emmanuel França, MD  
Brasil



Juliana Jordão, MD  
Brasil



Carlos R Antonio, MD  
Brasil



Ana P Meski, MD  
Brasil



Roseli Andrade, MD  
Brasil



Lincoln Fabricio, MD  
Brasil

# ETHEREA®: Versatility in Technology

**ATHENA**  
TECNOLOGIA PARA O TRATAMENTO ÍNTIMO FEMININO A LASER

**SPOT INLIFT**  
LASER LIFTING ATRAVÉS DE TRATAMENTO INTRAORAL

**IPL-SQ®**  
TECNOLOGIA IPL COM PULSO DIGITAL MICROPROCESSADO

**DUALMODE®**  
ER-YAG DE PULSO DUPLO, COM EFEITO COAGULATIVO COMBINADO

**INTENSE-IR®**  
TECNOLOGIA DE LUZ PARA SKIN TIGHTENING EM FACE E CORPO

**PRODEEP®**  
LASER FRACIONADO NÃO-ABLATIVO COM EFEITO SUBDÉRMICO

**ACROMA®**  
LASER Q-SWITCHED FRACIONADO E 2 COMPRIMENTOS DE ONDA

**LONGPULSE®**  
LASER ND:YAG DE PULSO VARIÁVEL - MICROPULSADO OU PULSO LONGO


**GOSMOOTH®**  
PADRÃO OURO EM RESURFACING FRACIONADO NÃO-ABLATIVO

# ETHEREA®: Versatility in Treatments

+ 70 TREATMENT INDICATIONS AVAILABLE AND GROWING

IPL-Sq®	intenseIR®	ACROMA-QS®	LongPulse®	ProDeep®	GoSmooth®	DualMode®
<ul style="list-style-type: none"> <li>• hair removal</li> <li>• pigmented and vascular lesions</li> <li>• anti-aging and wrinkles</li> <li>• active acne.</li> </ul>	<ul style="list-style-type: none"> <li>• <b>skin tightening on face and body areas</b></li> </ul>	<ul style="list-style-type: none"> <li>• LASER toning</li> <li>• tattoo removal</li> <li>• discrete pigmented and vascular lesions</li> </ul>	<ul style="list-style-type: none"> <li>• veins and vessels of the legs</li> <li>• hair removal in high phototypes</li> <li>• onychomycosis</li> <li>• LASER toning and thermal peeling for <i>flushing</i> ;</li> </ul>	<ul style="list-style-type: none"> <li>• Fractional non-ablative LASER</li> <li>• treatment of deep scars</li> <li>• stretch marks</li> <li>• melasma</li> <li>• hair therapy</li> </ul>	<ul style="list-style-type: none"> <li>• Fractional LASER for non-ablative FRACTIONAL RESURFACING</li> <li>• anti-aging and wrinkles</li> <li>• hair therapy</li> </ul>	<ul style="list-style-type: none"> <li>• Fractional ablative resurfacing</li> <li>• treatment of the intimate health of woman with LASER;</li> <li>• intimate rejuvenation</li> <li>• anti-aging and wrinkles</li> <li>• LASER peel for pigmentation on face, body and intimate region for woman</li> </ul>





# Handpieces



- 3, 5, 7 mm 1,064 nm
- 3.5 mm 532 nm
- 9/100 mtz/cm<sup>2</sup>


## INDICATIONS

- **removal of tattoos** with different colors and depths;
- **various pigmented lesions including**, post-inflammatory hyperpigmentation, ocher dermatites, hemosiderin deposition and keratoses;
- **removal of fine/light hair**
- **LASER *toning*** for all phototypes;
- the best therapeutic solution for **melasma**;

# LASER Nd:YAG 1,064 nm


*q-switched technology and photoacoustic effect*

ACROMA<sup>QS</sup>



**double rod design:** LASER rod generator and LASER rod amplifier: higher energy per shot

- **Nd:YAG 1064nm and KTP 532nm for multifunctional platform;**
- technology that seeks to add versatility and offer more treatments option;
- most LASERs use only **1 single rod for beam generation – less energy;**
- **ACROMA-QS<sup>®</sup> 2 rods, for generating and amplifying beam energy;**



## PARAMETERIZATION

- **energy:** 600 to 1500 mJ;
- **pulse duration** amount of energy delivered in a given time (20 ns);
- **frequency:** repetition rate, between shots;
- **spot size:** indicates which spot is in use. Larger spots sizes perform lighter treatments.

# ACROMA-QS®: Clinical Guide

## *LASER toning and rejuvenation*

ACROMA<sup>QS</sup>



Use of 1,064 QS Nd:YAG LASER, with collimated beam, low fluency and large spots;


- pore reduction;
- comedones;
- oiliness;
- texture and tonus;
- fine wrinkles;

- Very popular in Asia;
- multiple sessions, maximum safety and no downtime;
- No restriction to the phototype or tanning condition;

# ETHEREA-MX<sup>®</sup> ACROMA-QS<sup>®</sup>

*preventive maintenance and special care*

ACROMA<sup>QS</sup>



## Parameters for Black Peel:

1064nm 7mm spot

600mj

- care when cleaning especially **after the black peel procedure;**
- there can be no dirt on the lens because **this can irreversibly damage the cavity;**



# ETHEREA-MX<sup>®</sup> LASER GEMINI

*combined treatments*

ACROMA<sup>QS</sup> **New**

**LASER GEMINI:** combination of 2 wavelengths and identical half-actives, but with different pulse duration, for personalized effects;

**Nd:YAG q-switched, 1,064 nm**  
for skin spots and melanoses skin whitening  
in general  
**photomechanical effect only**

**Nd: YAG short-pulsed, 1,064 nm**  
for *toning, flushing* reduction, pore closure  
**thermal effect only**

- combination in the same treatment session;
- monthly or biweekly sessions;
- 6 to 8 treatment sessions for extended results;
- no *downtime*, restriction to phototype or tanning;

# ETHEREA-MX<sup>®</sup> LASER GEMINI

## *Treatment Protocol*

ACROMA<sup>QS</sup>

**New**

**LASER GEMINI: combination of 2 wavelengths in the same treatment session**

### **Acroma QS Handpiece**

**Nd:YAG q-switched, 1,064 nm**

**Spot 7mm , 1200mj, 3- 4 passes, until light to moderate Erythema**

### **LongPulse Handpiece**

**Nd: YAG short-pulsed, 1,064 nm**

**Dynamics Mode, Pulse 650 microsec**

**Tip 6mm, 10J/Cm<sup>2</sup>, 5.000 Pulses**

**Patient feedback**

**Frequency - up to 10HZ according to paciente tolerance**



# ETHEREA-MX<sup>®</sup> LASER GEMINI

*Courtesy Dr Valeria Campos*

ACROMA<sup>QS</sup>

**New**

**Before**



**After**



# ETHEREA-MX<sup>®</sup> LASER GEMINI

*Courtesy Dr Valeria Campos*

ACROMA<sup>QS</sup>

**New**

**Before**



**After**



# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>



- Long pulse ND:Yag 1064nm
- Gold standard for Vascular Lesions
- 4 spots – for improved results in different targets
- Spots plug and play
- LongPulse DYNAMICS, for *thermal peeling*



- 2, 3, 6 and 9 mm



- Efficacy of LASER Comes from a perfect combination of the parameters:

Spot size  
energy/fluency  
Pulse duration

## INDICATIONS

- **vascular lesions** in lower limbs (deep) and neck/face (superficial);
- **epilation**, especially for higher phototypes (Fitzpatrick V–VI)
- LongPulse DYNAMICS, for *thermal peeling* – **rejuvenation** (in general) and **onychomycosis**;

# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>

## *technical features*




<b>wave-length</b>	1,064 nm
<b>spot size</b>	2, 3, 6, 9 mm
<b>maximum fluence</b>	1,5 – 500 J/cm <sup>2</sup>
<b>pulse duration, LongPulse™</b>	3 to 60 ms
<b>pulse duration, Dynamics™</b>	300 μs to 1 ms

# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>

## *parametrization of the treatment*



QUESTION	ANSWER	ADJUSTMENT VARIABLE
WHAT IS THE TARGET, WHERE IS AND WHAT IS THE BARRIER TO BE TRANSPOSED BY THE LASER?	hemoglobin; HbO	wave-length
AT WHICH DEPTH?	very superficial, a few millimeters	wave-length; spot size
WHAT IS THE DIAMETER?	very thin; and 3 mm to 4 mm for reticular	pulse duration
WHAT IS THE ABSORPTION?	depends on the chosen wavel-ength	fluency



- **larger spot sizes** are indicated for the treatment of deeper lesions in **vessels and telangiectases of the lower limbs**;
- **smaller spot size** are indicated for the treatment of superficial lesions, **especially those originating in face**;

SPOT SIZE	PENETRATION
1 mm	0.8 mm
3 mm	1.5 mm
7 mm	3 mm
10 mm	4 mm
12 mm	4.5 mm
18 mm	5 mm

# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>

*parametrization of the treatment*



<b>LARGER DIAMETER VESSELS</b>	<b>LONGER PULSE DURATION</b>
<b>SMALLER DIAMETER VESSELS</b>	<b>SHORTER PULSE DURATION</b>

<b>STRUCTURE</b>	<b>THERMAL RELAXING TIME</b>
telangiectasias	10 ms at 20 ms
venulectasis	20 ms at 30 ms
reticular veins	30 ms at 60 ms

**PULSE DURATION OR TIME OF EXPOSURE TO ENERGY IS  
DETERMINED IN RESPECT OF THE VESSEL DIAMETER**



## Superficial Telangiectasias < 0,5mm

2mm 10ms 225-275J/cm<sup>2</sup>  
3mm 10-15ms 125-175J/cm<sup>2</sup>



## Telangiectasias from 0,5mm to 1mm

3mm 10-15ms 125-175J/cm<sup>2</sup>  
6mm 10-15ms 90-120J/cm<sup>2</sup>



# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>

*parametrization of the treatment*



## Venulectasis from 1 to 2mm

6mm 15-20ms 70-100J/cm<sup>2</sup>



## Reticular Veins from 2 to 3 mm

6mm 20-40ms 80-100J/cm<sup>2</sup>



# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>

*clinical reference guide: vascular lesions*



- **attention** to the phototype and tanned patients;
- remove creams and/or lotions before application;
- applications in areas with a lot of hair can epilate;
- observe clinical response and attention to the patient's pain;
- always initiate the tx by the reticularis before the telangiectasias;
- imperative the use of the epidermal cooler, SIBERIAN-FIT<sup>®</sup>;

# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>

*clinical reference guide: vascular lesions*



- **post-procedure:** apply local topical corticosteroids;
- **interval:** monthly, for the same area;

# Casebook and Results

## *vascular lesions*



**photo courtesy of:** Luiz Marcelo Viarengo. Vascular Surgeon. Jundiaía, SP. Brazil.

# Casebook and Results

## *vascular lesions*



treatment of telangiectases of the lower limbs with 1,064 LongPulse; 1 tx sessions, with intervals of 30 days

**photo courtesy of:** Rodrigo Kikuchi. Vascular Surgeon. São Paulo-SP. Brazil.

# Casebook and Results

## *vascular lesions*




treatment of telangiectases of the lower limbs with 1,064 LongPulse; 1 tx sessions, with intervals of 30 days

**photo courtesy of:** Rodrigo Kikuchi. Vascular Surgeon. São Paulo-SP. Brazil.

# Casebook and Results

*Vein Lakes: Courtesy Dr Rodrigo Kikuchi*



Single treatment  
Nd:YAG 1064nm

6mm  
40ms  
70J/cm<sup>2</sup>

The logo for etherea MX, featuring a colorful, multi-colored leaf-like icon above the text "etherea<sup>MX</sup>" in a lowercase, sans-serif font.

etherea<sup>MX</sup>

6mm 30-40ms 50-70J/cm<sup>2</sup>



# Casebook and Results

*Vein Lakes: Courtesy Dr Flavia Reginatto*




## Single treatment



6mm 20ms 90J/cm<sup>2</sup>


**Before**

**After**



# Casebook and Results

*Cherry: Courtesy Dr Rodrigo Kikuchi*



Single treatment  
Nd:YAG 1064nm

2mm  
10ms  
250J/cm



2mm 10ms 250-300J/cm<sup>2</sup>  
3mm 10-15ms 175-225J/cm<sup>2</sup>

# Casebook and Results


*hair removal*



**VYDENCE Medical Training Center**

# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>


*DYNAMICS MODE, micropulsed ND:Yag*



- LongPulse DYNAMICS, for *thermal peeling* – **rejuvenation** (in general); **onychomycosis, dilated pores, fine lines, scars, difused redness, rosácea flushing**
- **pulse duration**
  - 300 μs, 650 μs and 1 ms
- painless;
- no contraindication with respect to the phototype;

# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>

## *parameters and clinical guidelines*



- **1,000 to 2,000 shots** per quadrant, depending on the model;
- average of 6,000 to 10,000 shots on one whole face;
- if patient indicates very hot sensation , **reduce frequency of shots – 10, 7, 5, 2 Hz;**
- **immediate results**, but on average **1 to 8 treatment sessions** are performed ;
- **Interval:** 15 or 30 days;
- **Spot 6mm, 650micros, 10J/cm<sup>2</sup>**

# ETHEREA-MX<sup>®</sup> LongPulse<sup>®</sup>

*Dynamics Mode: Courtesy Dra Valeria Campos*



**Before**



**After 4 TX**



# ETHEREA-MX<sup>®</sup> DUALMODE<sup>®</sup> LASER Er:YAG 2,940 nm

DUALMODE



## INDICATIONS

- **Fractional Resurfacing;**
- **Photo Rejuvenation;**
- **Collagen Stimulation**, fine lines and moderate to severe wrinkles
- **Scars:** atrophic, hypertrophic, acne, traumatic and surgical;
- **Keratosis**
- **Benign pigmentary skin lesions**
- **Laser peel** with superficial tissue vaporization;
- **Warts.**

- Ø 2,5 and 6mm Collimated
- 8mm, 100 and 400Mtz/cm@

# ETHEREA-MX<sup>®</sup> DUALMODE<sup>®</sup> LASER Er:YAG 2,940 nm

DUALMODE



- **interval between sessions:** monthly or every 2 months;
- **average number of sessions:** 3 to 5;
- **special care** required in post-treatment

- Ø 2,5 and 6mm Collimated
- 8mm, 100 and 400Mtz/cm@



# ETHEREA-MX<sup>®</sup> DualMode<sup>®</sup>

## *technical features*

DUALMODE

<b>wave-length</b>	2,940nm
<b>spot size</b>	Collimated: 6mm Fractional: 8mm/100mtz and 8mm/400mtz
<b>maximum fluence</b>	Up to 60mj/mtz
<b>pulse duration</b>	300μs a 5 ms

# DualMode®: Clinical Guide

*fractional ablative resurfacing*

DUALMODE

<b>mJ/mtz energy</b>	the higher the greater the aggressiveness
<b>fluence mJ/cm<sup>2</sup></b>	the higher the greater the aggressiveness
<b>number of applications</b>	the higher the greater the aggressiveness
<b>pulse duration &lt;500 μs</b>	purely ablative effect
<b>pulse duration &lt;1ms</b>	coagulative effect
<b>double pulse, DualMode</b>	first ablative pulse followed by another coagulative


\* the longer the pulse duration, the greater the coagulative effect on the shot

\*\* DUALMODE pulse is more aggressive than the single shot

# LASER Er:YAG 2,940 nm

## *target-chromophore and absorption curve*

DUALMODE



- technology that revolutionized dermatology, introduced in 2004 by Mainstein et al\*
- works with the vaporization (ablation) of tissue
- there are 3 wave-lengths among fractional ablatives:
  - **2,790 nm:** Solid state LASER: YSGG (erbium-doped: yttrium-scandium-gallium-garnet)
  - **2,940 nm:** Solid state LASER, Er: YAG (erbium-doped: yttrium-aluminum-garnet)
  - **10,600 nm:** gas LASER, CO<sub>2</sub>

\* Manstein et al. FRACTIONAL PHOTOTHERMOLYSIS: A NEW CONCEPT FOR CUTANEOUS REMODELING USING MICROSCOPIC PATTERS OF THERMAL INJURY. LASERs Surg Med 2004;34:426-38.

# LASER Er:YAG 2,940 nm

*ablation, coagulation and ablation-coagulation*

DUALMODE


Ablative Skin Resurfacing  
(CO<sub>2</sub>, Er:YAG 2,940 nm)  
10–200 μm

Fractional Ablative  
Superficial Skin  
Resurfacing  
(CO<sub>2</sub>, Er:YAG 2,940 nm)  
10–70 μm

Fractional Non-  
ablative Skin  
Resurfacing  
600–1,000 μm

Fractional Ablative  
Skin Resurfacing  
100–300 μm

DualMode® Fractional  
Ablative/Coagulative  
Skin Resurfacing  
(CO<sub>2</sub>-like)  
100–300 μm




- Two types of Er:YAG – **LASER peel**, with pure ablative shots (*single-moded LASERS*), and **double pulse**, with combined ablative and coagulative shots (*dual-moded LASERS*);
- **DualMode® technology**: ablative and coagulative pulse in a single shot, simulating the coagulative effect inherent to CO<sub>2</sub> LASERS;

# LASER Er:YAG 2,940 nm

DUALMODE

*ablation, coagulation and ablation-coagulation*

Comparative	ABLATIVE	NON-ABLATIVE
<b>target-chromophore</b>	H2O	H2O
<b>chromophore absorption</b>	higher	lower
<b>direct effect</b>	vaporization	coagulation
<b>time-response</b>	48-72 hours.	24 hours.
<b>general advantages</b>	lasting results evident improvement high patient satisfaction less no. of sessions	satisfactory treatment extra-face tx safety reduced erythema downtime higher no. of indications greater versatility



## ETHEREA-MX<sup>®</sup> DualMode<sup>®</sup>: Results and Casebook

# LASER Er:YAG 2,940 nm

*After 2 Tx*

DUALMODE



# LASER Er:YAG 2,940 nm

*After 2 Txs*

DUALMODE





# LASER Er:YAG 2,940 nm

*After 1 Tx*

DUALMODE



# LASER Er:YAG 2,940 nm

*Courtesy: DR Valeria Campos – After 2 Txs*

DUALMODE



# LASER Er:YAG 2,940 nm

*Courtesy: DR Valeria Campos – After 2 Tx's*

DUALMODE



# ETHEREA-MX<sup>®</sup> GoSmooth<sup>®</sup>



- 8/100 mtz/cm<sup>2</sup>
- 10/400 mtz/cm<sup>2</sup>

- **Gold Standard Erbium-Glass 1540nm** technology of wide international acceptance in different global markets;
- **versatile:** *spots* available at 100 and 400 mtz/cm<sup>2</sup>;
- **better usability**, with square *spot*;
- optical arrangement system of the **high performance fractional lenses**;
- plus a fractional non-ablative LASER option for ETHEREA-MX;



- 8/100 mtz/cm<sup>2</sup>
- 10/400 mtz/cm<sup>2</sup>

## INDICATIONS

- fractional non-ablative *resurfacing*, stretch marks, scars – acne, surgical, recent, atrophic;
- **off-label:** melasma, alopecia and onychomycosis;
- better indication: **stretch marks** and **acne scars**;
- removal of superficial pigmented lesions – *skin toning*, whitening;

# ETHEREA-MX<sup>®</sup> ProDeep<sup>®</sup>



- **unique and exclusive** wave-length **in the market** - Nd-Yap 1,340 nm;
- **greater relative penetration** than the *gold standard* options available;
- coagulation effect at a more intense level, **greater bulking of subdermal heating**;

- 8/100 mtz/cm<sup>2</sup>
- 10/400 mtz/cm<sup>2</sup>
- 6 mm collimated



- 8/100 mtz/cm<sup>2</sup>
- 10/400 mtz/cm<sup>2</sup>
- 6 mm collimated

## INDICATIONS

- fractional non-ablative resurfacing, stretch marks, scars – acne, surgical, recent, atrophic;
- **off-label:** melasma, alopecia and onychomycosis;
- better indication: **stretch marks and acne scars;**
- **Clinical Studies:** Nodulocystic Acne resistant to isotretinoin. Hydrosadenite. Folliculitis Decalvant, Folliculitis Dissecans

# ETHEREA-MX<sup>®</sup> ProDeep<sup>®</sup>



- 8/100 mtz/cm<sup>2</sup>
- 10/400 mtz/cm<sup>2</sup>
- 6 mm collimated

- **interval between sessions:** monthly or every 2 months;
- **average number of sessions:** 3 to 7;
- Average downtime of **1 day with edema and slight erythema;**
- for **stretch marks**, relative to **3 days;**
- **Contraindication** to treatment: pregnancy;
- Frequently **combined** with **IPL Treatments.**



# ETHEREA-MX<sup>®</sup> ProDeep<sup>®</sup>


*technical features*




<b>wave-length</b>	Nd:YAP 1,340 nm
<b>spot size</b>	8 mm and 100 mtz/cm <sup>2</sup> 10 mm and 400 mtz/cm <sup>2</sup> 6 mm, collimated
<b>maximum fleunce</b>	200 mJ/mtz
<b>pulse duration</b>	3 to 20 ms
	with integrated cooling coupling; cold-air cooling

# ProDeep®: Clinical Guide

*How to treat*




Técnica Adequada



Técnica Não-Recomendada


- be careful with stacking during application;
- possible tissue retraction if the mode of application is not followed;
- subdermal effect may generate unwanted heating bulking



Técnica Não-Recomendada

# LASER Nd:YAP 1,340 nm

*ablation, coagulation and ablation-coagulation*



- ...
- Flg 2-C: shows depth (612,8 µm) and width/diameter (418,0 µm) of fractional microscopic treatment zone at 1,340 nm (fractional non-ablative LASER) in the tissue at D3 with 8 x 8 mm spot size, 100 MTZ/cm<sup>2</sup> and 200 mJ/MTZ.

# Casebook e Resultados

## *rejuvenescimento*



foto cortesia de Valéria Campos, MD. Jundiaí, SP, Brasil.

Workshop ETHEREA 2014

fotos de antes e depois de 1 sessão de tratamento

100 mtz/cm<sup>2</sup>, 90mJ/mtz, 5 ms;

# Casebook e Resultados

*acne*



foto cortesia de Valéria Campos, MD. Jundiaí, SP, Brasil.

Workshop ETHEREA 2014

fotos de antes e depois de 1 sessão de tratamento

100 mtz/cm<sup>2</sup>, 120mJ/mtz, 5 ms

Confidencial. Uso interno. Distribuição e/ou reprodução proibidos.




# “Stretch Marks treatment with non-ablative technologies”

**Dra. Luiza Pitassi**

Doutorado e Mestrado - Dermatologia - UNICAMP


Coordenadora do Ambulatório de Cosmiatria - UNICAMP



Faculdade de Ciências Médicas  
Dermatologia Unicamp



# Non-ablative Lasers



*Foto: arquivo pessoal Dra. Luiza Pitassi*

Tecnologias NÃO ABLATIVAS	Comprimentos de onda
Er:Glass	1540 nm
Nd:YAP	1340 nm
IPL	400-1200nm

1-Taub AF. Fractionated delivery systems for difficult to treat clinical applications: acne scarring, melasma, atrophic scarring, striae distensae, and deep rhytides. *J Drugs Dermatol.* 2007; 6(11):1120-8.  
2-El Ammar, A. B. P. C., Campos, V. B., & Duquia, R. P. (2010). Análise clínica de 16 pacientes consecutivos tratados com LUX 1540® para rejuvenescimento. *Surgical & Cosmetic Dermatology*, (3), 221-224.

# ∞ Non-ablative Lasers

The target is **Water**

**SAFE FOR ALL SKIN TYPES!!!**



Clark CM, Silverberg JI, Alexis AF . A retrospective chart review to assess the safety of nonablative fractional laser resurfacing in Fitzpatrick skin types IV to VI. J Drugs Dermatol. 2013 Apr;12(4):428-31.



## 1540 nm



8mm - 100mtz/ cm<sup>2</sup>

- Min Fluence- 15 mJ
- Max Fluence– 95mJ
- **55-65-75mJ/ 10ms**

## 1340 nm


8mm - 100mtz/ cm<sup>2</sup>

- Min Fluence- 60 mJ
- Max Fluence– 200 mJ
- **110mJ/ 3ms**




# LASER Nd: AP 1,340 nm

## target-chromophore and absorption curve




- there are several wave-lengths among the fractional non-ablative;
- gold standard technologies for rejuvenation (commercial) are the main ones today;
  - **1,320 nm:** Solid state LASER, Nd: YAG;
  - **1,440 nm:** Diode LASER, by fiber;
  - **1,340 nm:** Solid state LASER, Nd:YAP;
  - **1,540 nm:** Solid state LASER, Er:Glass;
  - **1,540 nm:** Solid-state laser, Er:Glass, with optical fiber delivery and scanner;

# Nd YAP:1340nm



**1340 – 8mm / 60J**



**1340 – 8mm / 130J**

*Fotos: arquivo pessoal Dra. Luiza Pitassi*

# Technologies Association

## IPL + LASER

### RED STRECH MARKS

- IPL + Fractional Non-Ablative Laser + optimal results
- 3 to 5 treatments 30-60 day apart

➤ ***IPL: 540nm, 13- 16mJ/15ms, 1 pass (Skin Thypes I-IV)***

➤ ***ProDeep 1340nm: 100-140mJ/3ms***

***Spot 8mm/100 mtz /cm2***

# ETHEREA-MX<sup>®</sup> ProDeep<sup>®</sup>

Results: Courtesy of Dr Carlos Roberto Antonio

## Nodulocystic Acne resistant to isotretinoin

### Protocol and Guidelines:

2 to 6 Tx's

30 days apart

Spot 8mm/100mtz

100mj / 3ms



FIGURA 1: Lesões acneicas antes do tratamento (A) e após o tratamento (B)



FIGURA 2: Lesões acneicas antes do tratamento (A) e após o tratamento (B)



FIGURA 3: Lesões acneicas antes do tratamento (A) e após o tratamento (B)

Hidradenitis : Courtesy of Dr Carlos Roberto Antonio

## Protocol and Guidelines:

4 Txs, 30 days apart

Spot 8mm/100mtz

100mj / 3ms

**Before**



**After 4 months**



**3 years FUP**



# CICATRIZ DE ACNE

PRODEEP

Acne Scars: after 3 Txs



Antes



Depois

PONTEIRA: 100 mtz/cm<sup>2</sup>  
ENERGIA: 125 mJ/mtz

TEMPO DE PULSO: 5 ms  
NUMERO DE PASSES: 1



Uma empresa  
Adavium Medical

# CICATRIZ DE ACNE

PRODEEP



**Antes**



**Depois**

PONTEIRA: 100 mtz/cm<sup>2</sup>  
ENERGIA: 125 mJ/mtz

TEMPO DE PULSO: 5 ms  
NUMERO DE PASSES: 1

 **vyidence**  
medical

Uma empresa  
Adavium Medical



# CICATRIZ DE ACNE

PRODEEP



**Antes**



**Depois**

PONTEIRA: 100 mtz/cm<sup>2</sup>  
ENERGIA: 120 mJ/mtz

TEMPO DE PULSO: 5 ms  
NUMERO DE PASSES: 1

 **vyidence**  
medical

Uma empresa  
Adavium Medical

# CICATRIZ DE ACNE

PRODEEP



**Antes**



**Depois**

PONTEIRA: 100 mtz/cm<sup>2</sup>  
ENERGIA: 125 mJ/mtz

TEMPO DE PULSO: 5 ms  
NUMERO DE PASSES: 1

# CICATRIZ DE ACNE

PRODEEP



**Antes**



**Depois**

PONTEIRA: 100 mtz/cm<sup>2</sup>  
ENERGIA: 125 mJ/mtz

TEMPO DE PULSO: 5 ms  
NUMERO DE PASSES: 1

# LASER Nd: AP 1,340 nm

*Scars: Courtesy Dr. Bruna Lenz*




**Before**



**4 Txs - 30 days apart**  
**Spot 8/100mtz**  
**115- 140mj/mtz**  
**5ms**

**After**





- **highly flexible system:** can treat a set of conditions with a single device;
- **effective:** the results are good, similar to some laser applications;
- **safe:** proven by hundreds of scientific papers and years of clinical practice;



- *all-in-one system*: **one handpiece, with multiple plug-and-play cutting filters**, slice, at 400, 540, 580, 640 and 695 nm;
- exclusive 695 nm cutting filter, **indicated for high phototypes**;
- pulses of 5 to 100 ms, which offer **aggressiveness when needed and safety when required**;
- **a single handpiece that offers unparalleled versatility** in the market, maximizing the potential indications of technology;




- **wide coverage area**, with 40 x 12 mm, 4.8 cm<sup>2</sup>;
- VascuTips<sup>®</sup> with Ø 8 mm and 12 x 12 mm **for localized lesions**;
- repetition rate of up to 2 Hz, offering **much faster treatment**;
- **wide coverage area means greater profitability**, with fewer shots needed to treat the same area;



- **sapphire contact cooling**, ranging from 5 to 20 °C with 5 levels of intensity;
- sapphire is different from other materials because of its **high efficiency cooling**;
- protects the epidermis against possible collateral damage to the treatment;
- mild anesthetic effect, reducing patient discomfort during the treatment;
- allows higher fluences to be used for greater and better results;





- fluency up to 33 J/cm<sup>2</sup>, with delivery of Square-Wave Pulse wave, optimizing safety and increasing results;
- on the free discharge wave, **the risk of adverse effects is higher**, since the peak temperature reached is much higher than that actually required;
- **in the format wave Square-Wave Pulse, the values are constant, with calculated and controlled discharge**, ensuring uniform distribution of energy throughout the pulse;

### *Pigmented lesions*

- **540 nm:** light phototypes, clearer and/or superficial lesions;
- **580 nm:** high phototypes, very dark and/or deep lesions;
- **IMPORTANT:** when using VascuTips, maintain level of cooling in 5
- Late reactions to treatment, considered desirable are moderate erythema around the lesion, with corresponding darkening, with no visible change in the adjacent tissue (localized effect on the lesion);

### *Vascular lesions*

- **540 nm:** light phototypes, superficial lesions;
- **580 nm:** high phototypes, deep lesions;
- **IMPORTANT:** fluency and pulse time should be adjusted according to vessel size and patient's phototype;
- when using vascular tips, maintain cooling at level 5.
- do not press the vessels before and during the shot, since this implies the blood dispersion and consequent treatment ineffectiveness;
- do not treat vascular lesions superimposed by pigmented lesions and/or hair;
- indications for face and trunk only, due to the superficial location of the lesions in the tissue (lower limbs lesions tend to deepen);

### Active Acne

- **filter 640 nm:** anti-inflammatory activity
- **filter 400 nm:** activity in protoporphyrin produced by *P. Acnes*;

<b>level of cooling</b>	levels 4 or 5 for comfort and safety
<b>number of sessions</b>	from 1 to 10 *depends on degree of involvement and lesions
<b>interval between sessions</b>	7 to 14 days * do not exceed more than 14 days
<b>time for results</b>	the expected improvement is usually immediate
<b>relative downtime</b>	none

### *Rejuvenation*

- for cases of rejuvenation associated with **pigmented lesions** and/or **vascular lesions**, use parameters corresponding to the indication and only after that proceed with the rejuvenation treatment;
- the **640 and 695 filters** can be used in combination with all other treatment modalities, offering an interesting complementary alternative to the procedures;



Relato de Caso

## IV Encontro dos Usuários do ETHEREA



Uma empresa  
Adavium Medical

# Mancha Vinho do Porto

Luz Pulsada – filtro 540 – ETHEREA



Antes



Depois

*J.M.J.S.M, 37 anos, com diagnóstico de Mancha Vinho do Porto, iniciou tratamento com Luz Pulsada (filtro 540-Etherea), foto após 5 sessões, realizadas com intervalo mensal, com parâmetros que variaram 17-19 J/cm e 15 ms. Nota-se clareamento importante, além da melhora das lesões nodulares. Fotos cortesia de: Dra Heloisa B Tiraboschi Romanholi*

## Active Acne: Courtesy Dra Celia Kalil



Before

After 10 Txs

640nm 16J/cm<sup>2</sup> 30ms Face

640nm 18J/cm<sup>2</sup> 100ms Cheeks and chin

400nm 12J/cm<sup>2</sup> 30ms Cheeks and chin

## Poiquilodermia de Civatte: Courtesy Dra Celia Kalil



**Before**

**After 2 Txs**

**IPL 540nm  
15ms 13 - 15J**



# IPL + Prodeep 1340nm: Courtesy Dra Célia Kalil



30/30d

540nm 12x12mm 10ms 15J - 20/cm<sup>2</sup>

Prodeep 1340nm 8mm 100MTZ/cm<sup>2</sup> 5ms 70-90mJ/MTZ

# Casebook e Resultados

## *lesões pigmentares*



fotos de antes e depois do tratamento

Confidencial. Uso interno. Distribuição e/ou reprodução proibidos.

# Casebook e Resultados

## *lesões pigmentares*



fotos de antes e depois do tratamento

# Casebook e Resultados

## *lesões vasculares*

ipl·sq



fotos de antes e depois do tratamento

Confidencial. Uso interno. Distribuição e/ou reprodução proibidos.

# Casebook e Resultados

## *lesões vasculares*

ipl·sq



fotos de antes e depois do tratamento

# Casebook e Resultados

## *cicatrizes*

IPL-Sq



fotos de antes e depois de 5 sessões de tratamento

- 1: 540 nm 16 J/cm<sup>2</sup> 20 ms [IPL-Sq]; 100 mtz/cm<sup>2</sup>, 130mj/mtz, 3 mm [ProDeep];
- 2: 540 nm 17/18 J/cm<sup>2</sup> 15 ms; [IPL-Sq]; 100 mtz/cm<sup>2</sup>, 140mj/mtz, 3 mm [ProDeep];
- 3-5: 540 nm 17/18 J/cm<sup>2</sup> 10 ms; [IPL-Sq]; 100 mtz/cm<sup>2</sup>, 150mj/mtz, 3 mm [ProDeep]

# Casebook e Resultados

## *lesões vasculares*



foto cortesia de Lilian Braga, MD. Sorocaba, SP, Brasil

fotos de antes e depois de 5 sessões de tratamento  
540 nm 15 J/cm<sup>2</sup> 20 ms;

# Casebook e Resultados

*lesões vasculares + lesões pigmentares*

ipl·sq



fotos de antes e depois do tratamento



# Casebook e Resultados

## *cicatrizes*



**foto cortesia** de Célia Kalil, MD. Porto Alegre, RS, Brasil.


IV Encontro dos Usuários ETHEREA

fotos de antes e depois de 4 sessões de tratamento

540 nm 15 J/cm<sup>2</sup> 15 ms;

# Casebook e Resultados

## *rejuvenescimento*



**foto cortesia** de Airá Novello Villar, MD. Concórdia, SC, Brasil.

III Encontro dos Usuários ETHEREA

fotos de antes e depois do tratamento

540 nm 15 J/cm<sup>2</sup> 15 ms; [IPL-Sq]

160mJ/mtz 5 ms; [ProDeep]

Confidencial. Uso interno. Distribuição e/ou reprodução proibidos.

# ETHEREA-MX<sup>®</sup> Intense IR<sup>®</sup>

*wavelength and absorption curve*

INTENSE IR




- **IntenseIR** source of high energy infrared light
- Gold Standart for **Skin Tightening**;
- **Face, Neck** and **body** treatments;
- **no contraindication** with respect to the **phototype**;
- Painless (dynamics mode);
- **No downtime.**

# IPL: Intense Pulsed Light

## *wavelength and absorption curve*

INTENSE IR



- Thermal image shows the heating in the target tissue. Adjacent zones are protected and cooled. The epidermis is maintained at a safe temperature (less than 40 ° C).

# ETHEREA-MX<sup>®</sup> Intense IR<sup>®</sup>

*interface parameterization and usability*

INTENSE IR



- **fluence:** in J/cm<sup>2</sup>, referring to energy delivered in the area of the tx;
- **pulse time:** in seconds, time that the tissue will be exposed to IR.
- **Pre-cooling:** in seconds, before the pulse. Protects the epidermis, preventing subsequent damage(only stamp mode)
- **Post-cooling:** in seconds. After the pulse(only stamp mode)

# ETHEREA-MX<sup>®</sup> Intense IR<sup>®</sup>

## *Treatment Protocol*

INTENSE IR



- **3 treatments, 30 days apart. 6 months maintenance.**
- **Dynamics application.**
- **5 – 8 pulses per area 20 to 30 on face, and 25 to 35 on the neck.**

# ETHEREA-MX<sup>®</sup> Intense IR<sup>®</sup>

## Casebook: Skin Tightening

INTENSE IR



TP 10 s, 91 J/cm<sup>2</sup>, pre/post-cooling 0

# ETHEREA-MX<sup>®</sup> Intense IR<sup>®</sup>

## Casebook: Skin Tightening

INTENSE IR



TP 10 s, 91 J/cm<sup>2</sup>, pre/post-cooling 0

Confidential and proprietary



# ETHEREA-MX<sup>®</sup> Intense IR<sup>®</sup>

## Casebook: Skin Tightening

INTENSE IR



TP 10 s, 91 J/cm<sup>2</sup>, pre/post-cooling 0

Confidential and proprietary

# ETHEREA-MX<sup>®</sup> Intense IR<sup>®</sup>

## Casebook: Skin Tightening

INTENSE IR



TP 10 s, 110,5 J/cm<sup>2</sup>, pre/post-cooling 0, cooling nível 5

RELEVANT PLATFORM FEATURES	MX	JOULE	XEO	ICON	M22	HRMNY
IPL all-in-one handpiece	✓	✓	✗	✗	✓	✗
IPL square-wave pulse	✓	✓	✓	✓	✗	✗
IPL accurate spot sizes	✓	✓	✗	✓	✓	✗
halogen IR body & face tx options	✓	✗	✓	✓	✗	✓
q-switched LASER handpiece	✓	✗	✗	✗	✓	✓
QS fractional q-switched system	✓	✗	✗	✗	✗	✓
1064 higher energy (vascular)	✓	✓	✓	✓	✓	✓
1064 therma-peel, short-pulsed	✓	✓	✓	✗	✗	✗
1340 unique, deep epidermal	✓	✗	✗	✗	✗	✗
1540 gold standard LASER	✓	✗	✗	✓	✓	✓
2940 dual-effect Er:YAG (CO <sub>2</sub> -like)	✓	✓	✗	✓	✗	✗
2940 short pulsed LASER peel	✓	✓	✗	✓	✗	✗
2940 intraoral LASER lift	✓	✓	✗	✗	✗	✗
2940 intimate LASER tx option	✓	✓	✗	✗	✗	✗



VYDENCE® Medical Office US 3340 Hillview

Avenue Palo Alto, CA. United States

VYDENCE® Medical Facility R&D Rua Aldo  
Germano Klein, 359. CEAT

13573-470 São Carlos, SP. Brazil

VYDENCE® Medical Training Center

Av Fagundes Filho, 486, SI 96. S Judas  
04304-000. São Paulo-SP. Brazil

ADAVIUM® Medical Office BRA Av.

Queiroz Filho, 1560. T Gaivota  
05319-000. São Paulo-SP. Brazil



Uma empresa  
Adavium Medical